

Product Name :
Blood Bag Sealer

Product Code :
NLABMEDICAL-M3200008

Description :

Blood Bag Sealer India, Blood Bag Sealer manufacturer India, suppliers, Medical Blood Bag Sealer India, Hospital Blood Bag Sealer India

Technical Specification :

Sealing operation is activated automatically when the tube is based in between the electrodes.

Electrodes are made from high grade metal for perfect sealing. Caps provided in front of electrodes are removable, so it is easy to clean.

LED's are provided for indications of Power, Read, Seal, Ground.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003