

Product Name :
Blood Bag Tube Sealer

Product Code :
NLABMEDICAL-M3200002

Description :

Blood Bag Tube Sealer India, Blood Bag Tube Sealer manufacturer India, suppliers, Medical Blood Bag Tube Sealer India, Hospital Blood Bag Tube Sealer India

Technical Specification :

Oscillation Frequency: 40.68 Mhz

Radio Frequency output: 15W

Sealing time: Approx. 1.5 sec.

Input Volt: 230V AC

Power Frequency: 50 Hz

Power Consumption: 100 VA

Weight: 5.5 Kg approx

Size: 180 mm x 215 mm x 245 mm (HxDxW)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003