

<b>Product Name :</b> Blood Bank Freezers	<b>Product Code :</b> NLABMEDICAL-M3200009
<b>Description :</b> Blood Bank Freezers India, Blood Bank Freezers manufacturer India, suppliers, Medical Blood Bank Freezers India, Hospital Blood Bank Freezers India	
<b>Technical Specification :</b> Temperature control  Microprocessor control and electrical boosting system ensuring a constant temperature rang $4^{\circ}\text{f} \pm 1^{\circ}\text{f}$ .  Digital display of average temperature upper and lower set points with resolution and monitoring of cabinet temperature  Adjustable temperature control range $+2^{\circ}\text{f}$ to $+8^{\circ}\text{f}$ and increment $0.5^{\circ}\text{f}$ , factory pre-set to $+4^{\circ}\text{f}$ .  Safety control system  Dual display of operation information ( digital display and recorder)  Alarm for five defects (over temperature and under temperature, power failure, sensor defect, door opening, low voltage of battery backup)  Three warning performances ( audible/visual and remote signal)  Dual system of temperature data recording when power fails (digital display and recorder)  Dual safety operation programs (sensor defect and disorganized digits)  Ergonomic design  Safety lock to prevent accidental access  Large LED for temperature information  Sorted blood storage to ease accesses  Wide voltage for 187~242V power supply.	

Wide climatization for operation at 10â„ƒ~35â„ƒ ambient temperature.

Transparency design

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003