

**Product Name :**  
Blood Bank Freezers

**Product Code :**  
NLABMEDICAL-M3200009

**Description :**

Blood Bank Freezers India, Blood Bank Freezers manufacturer India, suppliers, Medical Blood Bank Freezers India, Hospital Blood Bank Freezers India

**Technical Specification :**

Temperature control

Microprocessor control and electrical boosting system ensuring a constant temperature range  $4\pm 1$ .

Digital display of average temperature upper and lower set points with resolution and monitoring of cabinet temperature

Adjustable temperature control range  $+2$  to  $+8$  and increment  $0.5$ , factory pre-set to  $+4$ .

Safety control system

Dual display of operation information ( digital display and recorder)

Alarm for five defects (over temperature and under temperature, power failure, sensor defect, door opening, low voltage of battery backup)

Three warning performances ( audible/visual and remote signal)

Dual system of temperature data recording when power fails (digital display and recorder)

Dual safety operation programs (sensor defect and disorganized digits)

Ergonomic design

Safety lock to prevent accidental access

Large LED for temperature information

Sorted blood storage to ease accesses

Wide voltage for 187~242V power supply.

Wide climatization for operation at 10?~35? ambient temperature.

Transparency design

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003