

**Product Name :**  
Blood Bank Refrigerator Deluxe Model

**Product Code :**  
NLABMEDICAL-M3200011

**Description :**

Blood Bank Refrigerator Deluxe Model India, Blood Bank Refrigerator Deluxe Model manufacturer India, suppliers, Medical Blood Bank Refrigerator Deluxe Model India, Hospital Blood Bank Refrigerator Deluxe Model India

**Technical Specification :**

Blood Bank Refrigerators are available in various capacities, designed to fit a variety of space and storage requirements.

The refrigerators are equipped with heavy-duty compressors to maintain a vital operation temperature of 4.0°C Superior refrigeration technology and precise temperature control to ensures reliable blood storage.

Inner and outer chamber both made of Stainless Steel (SS-304). Door, panel, base plate and top cover made of C.R.C. duly powder coated.

Stainless Steel Baskets are provided for storage inside.

Foamed - in - place CFC Free PUF insulation ensures temperature stability and reduced energy temperature.

Stainless steel heavy duty basket ensures even temperature distribution.

Micro - controller controls the internal temperature between 4°C to 6°C with an accuracy of  $\pm 0.5^{\circ}\text{C}$ , for storing blood bottles or blood bags.

The cooling is effected by a hermetically sealed CFC Free Eco Friendly Compressor.

To work on 220/230 volts A.C

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

