

**Product Name :**  
Blood Coagulator Analyzer

**Product Code :**  
NLABMEDICAL-M3200012

**Description :**

Blood Coagulator Analyzer India, Blood Coagulator Analyzer manufacturer India, suppliers, Medical Blood Coagulator Analyzer India, Hospital Blood Coagulator Analyzer India

**Technical Specification :**

Detector channels: 1

Position for Reagents: 2

Position for Sample: 5

Light source: 470nm long-life cold light source

Sample consumption: 20-40uL

Mini-reagent consumption: 20uL

Homoeothermic controlling: 37 $\pm$ 0.3 $^{\circ}$ C

Repeatability: CV $\leq$ 2%

Linear error: r $\geq$ 0.98

Measurement startup: online adding, auto timing, manual startup system

Memory: 10000 results

Counting function: can be converted on S, %, PTR, INR, mg/dL

Printer: built in thermal printer

Display: 192\*64 LCD display

Size: 280mm(L)\*300mm(B)\*160mm(H)

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003