

Product Name :
Blood Donor Chair

Product Code :
NLABMEDICAL-M2700001



Description :

Blood Donor Chair

Technical Specification :

The frame is made of Stainless Steel Round Tube.

Leg, Seat, Back and Armrest padded and covered with washable plastic material.

Simultaneous movement of back and leg rest.

Adjustable armrest.

Foldable Size : 100 L x 55 W x 120 H cms

Open Size 175 L x 55 W x 50 H cms

We are Manufacturer, Exporters and Suppliers of Blood Donor Chair at tender and export competitive prices from India.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments