

**Product Name :**  
Body Basic CPR Manikin Advanced Along With Printer

**Product Code :**  
NLAB-MEDICALMODELAB11031



**Description :**

Body Basic CPR Manikin Advanced Along With Printer

**Technical Specification :**

This model is designed according to the international CPR standard.

In this model is inbuilt light indicator on the printer for the correct & incorrect.

Artificial respiration (mouth-to mouth resuscitation & incorrect compression on the chest.

This model is designed according to the international CPR standard.

In this model is inbuilt light indicator on the printer for the correct & incorrect.

Artificial respiration (mouth-to mouth resuscitation & incorrect compression on the chest.

The inbuilt LCD on the printer shows the electrocardiogram for the experiment done.

The model also shows whether the operation is correct or not by light indicator, signal indicator, language indicator and grade printing.

In this when the patient revives the cardiogram will appear.

Features:-

High durability

Light weight

Easy movable

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003