

Product Name :
Body Basic CPR Manikin Female with Monitor and Printer

Product Code :
NLAB-MEDICALMODELAB11026



Description :

Body Basic CPR Manikin Female with Monitor and Printer

Technical Specification :

This model is designed according to the international CPR standard.

Simulate the patient revivals through first aid by measures of heart press on the chest and mouth respiration.

As soon as the patient revivals, his neck artery and heart recover pulsate continuously, his enlarged pupils shrink to normal.

The model also shows the operation is correct or not by singal indication.

Language indication and garde printing.

An ideal model for use in hospitals, medical colleges and schools.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

