

Product Name :
Bowl Stand Double Metal Base

Product Code :
NLABMEDICAL-M2810004



Description :

Bowl Stand Double Metal Base

Technical Specification :

Bowl Stand Double Metal Base

The pipe is made of Metal.

It comes with Base having 5 cms castors.

The base is made of Metal.

It comes with a double bowl.

The bowl is made of stainless steel and the size is 35 cms.

Knock Down Construction to save freight.

We are manufacturers, Exporters, and Suppliers of Bowl Stand Double Metal Base at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003