

**Product Name :**  
Bowl Stand Double Metal Base

**Product Code :**  
NLABMEDICAL-M2530007



**Description :**

Bowl Stand Double Metal Base

**Technical Specification :**

Specification :-

Pipe is made of Metal.

It comes with Base having 5 cms castors.

Base is made of Metal.

It comes with double bowl.

The bowl is made of stainless steel and the size is 35 cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Bowl Stand Double Metal Base at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003