

Product Name :
Bowl Stand Double Plastic Base

Product Code :
NLABMEDICAL-M2530006



Description :

Bowl Stand Double Plastic Base

Technical Specification :

Specification :-

Pipe is made of Metal.

It comes with Base having 5 cms castors.

Base is made of PU (Polyurethane).

It comes with double bowl.

The bowl is made of stainless steel and the size is 35 cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Bowl Stand Double Plastic Base at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003