

Product Name :
BP Machine Digital Deluxe

Product Code :
NLABMEDICAL-M2480001



Description :

BP Machine Digital Deluxe

Technical Specification :

BP Machine Digital Deluxe

Display : LCD Digital Display

Display Size : L 70mm x W 60mm

Display Backlight : Yes Orange Color

Body Color : White and Orange

Measurement Method : Oscillometric Method

Measurement Range

Pressure : 0 to 280 mmHg

Pulse : 40 to 180/min

Accuracy / Calibration

Pressure : ± 3 mmHg or 2% of reading

Pulse : +5% of reading

Inflation : Automatic by Electric Pump

Deflation : Automatic Pressure release valve

Rapid Air Release : Automatic exhaust valve

Power Source : 4 "AA" Batteries 1.5V, USB Cable

Battery Life : Approx. 700 measurements (using new

Alkaline Batteries) Automatic Power-off Function

Memory : Upto 90 measurements with date and time for each

Arm Cuff Size : 22-36 cm Medium

Functions:

WHO Indicator, Arrythmia Checking, Average Function

Date, Time and Room Temperature Display when not in use.

Device Dimension : W 105mm x L 130mm x H 60mm

Contents with Product

1 x Bp Box

1 x Bp Monitor

1 x Bp Cuff 22-36 cm

1 x Product Manual

4 x "AA" Battery

1 x USB Cable

1 x Storage Pouch

We are Exporters and Suppliers of Desco BP Machine Digital Deluxe at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003