

**Product Name :**  
Brazilian Test Apparatus

**Product Code :**  
NLAB-COMPSLAB280003



**Description :**

Brazilian Test Apparatus

**Technical Specification :**

Diameter: 50 mm - 100 mm

Automation Grade: Automatic

Pressure gauge capacity: 0-200 KN x 2 KN

Material: Stainless Steel

Display: Analog

The specimen is held in circular jaws, this is primarily similar to a compression Machine and consists of a small load frame having a sturdy base with two vertical threaded rods and an adjustable crosshead. The hydraulic jack is fitted at the center of the base of the load frame. The jack of the load frame is self-retracting and two plain platens are supplied.

A maximum pointer is also provided on the gauge a pair of semi-circular platens for 50mm dia, samples also provided. This instrument can be used in the field also.

Optional Extras: Pairs of jaws for sample dia.

- (a) 60mm
- (b) 70mm
- (c) 80mm
- (d) 90mm
- (e) 100mm.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003