

Product Name :
Brinell Hardness Tester

Product Code :
NLAB-MECHANICALAB400001



Description :

Brinell Hardness Tester

Technical Specification :

Analog readout ;

Height of tests: 20 mm (minimum);

Depth of throat : 17 mm (minimum);

Preload 10 kg;

Comes With:

Brinell indenter 2.5 mm ball

Protective cover

Accessories:

2 diamond penetrators

penetrators 2 1/16 " ball

2 • 2.5 mm ball indenter

1 standard Rockwell

1 anvil plate

1 anvil wide V

V 1 anvil in close

one set of ball 1/16 " recharging

1 set of ball Brinell

Related Services:

Switching on and testing.

Providing necessary for the maintenance of the machine

Detailed English language manual and maintenance.

User training for one (01) day minimum.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003