

**Product Name :**  
Brownlee Electrolysis Apparatus

**Product Code :**  
NLAB-EDUSCHOOLEQ143001



**Description :**

Brownlee Electrolysis Apparatus

**Technical Specification :**

Brownlee Electrolysis Apparatus

The most popular form of Electrolysis Apparatus.

Consists of Electrodes attached to insulated connecting rods which are in-turn connected to terminals.

Mounted on a non-conducting plastic holder, designed to rest on top of a glass jar / beaker.

The plastic holder has clips to hold two test tubes.

Operates on 6 - 10v DC, power supply not included.

**Supplied complete with:-**

Beaker (1000mL)

Graduated Test Tubes

Un-graduated Test Tubes

Rubber Stoppers for Test Tubes

A preview of the included instructions can be found at the bottom of this listing.

Measurements are approximate.

Instructions Included.

Shipping Weight 5.00 lbs

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003