

**Product Name :**  
BT Bottle Top Dispenser

**Product Code :**  
NLAB-ENGINEERINGLB97002



**Description :**

BT Bottle Top Dispenser

**Technical Specification :**

BT Bottle Top Dispenser

Six unique volume ranges:

0.25 - 2.5 ml

b. 0.5 - 5 ml

c. 1 - 10 ml

d. 2.5 - 30 ml

e. 5 - 60 ml

f. 10 - 100 ml

Fully Autoclavable at 121°C, 15psi. for 10 - 15 mins

Re-Circulation Valve

Prevents loss of reagent during purging.

#### Adjustable Delivery Nozzle

Adjustable delivery nozzle to facilitate dispensing ease in all demanding laboratory conditions.

#### Specially Designed Volume Adjustment Knob

180° rotation for easy & effortless volume setting.

#### 360° Rotation

Specially designed adaptor allows 360° rotation, providing full flexibility in working.

#### Springless Valve

Specially designed, springless, PTFE valve manifold ensures high chemical compatibility and jamming free valve functioning.

#### Unique Piston

PTFE piston with ETP O-Ring ensures very high chemical compatibility & wiper like smooth, effortless piston movement.

#### Calibration

Specially designed calibration tool is provided for convenient and quick in-lab user re-calibration. This is in compliance with GLO/ISO norms.

#### Adaptors

Adaptors to fit most of the laboratory reagent bottles are provided in the following sizes: 28, 30,32, 36, 40, and 45 mm.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003