

Product Name :
BT Bottle Top Dispenser

Product Code :
NLAB-ENGINEERINGLB97002



Description :

BT Bottle Top Dispenser

Technical Specification :

BT Bottle Top Dispenser

Six unique volume ranges:

0.25 - 2.5 ml

b. 0.5 - 5 ml

c. 1 - 10 ml

d. 2.5 - 30 ml

e. 5 - 60 ml

f. 10 - 100 ml

Fully Autoclavable at 121°C, 15psi. for 10 - 15 mins

Re-Circulation Valve

Prevents loss of reagent during purging.

Adjustable Delivery Nozzle

Adjustable delivery nozzle to facilitate dispensing ease in all demanding laboratory conditions.

Specially Designed Volume Adjustment Knob

180° rotation for easy & effortless volume setting.

360° Rotation

Specially designed adaptor allows 360° rotation, providing full flexibility in working.

Springless Valve

Specially designed, springless, PTFE valve manifold ensures high chemical compatibility and jamming free valve functioning.

Unique Piston

PTFE piston with ETP O-Ring ensures very high chemical compatibility & wiper like smooth, effortless piston movement.

Calibration

Specially designed calibration tool is provided for convenient and quick in-lab user re-calibration. This is in compliance with GLO/ISO norms.

Adaptors

Adaptors to fit most of the laboratory reagent bottles are provided in the following sizes: 28, 30,32, 36, 40, and 45 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003