

Product Name :
Buchner Funnel Porcelain Deluxe

Product Code :
NLAB-EDUSCHOOLEQ166002



Description :

Buchner Funnel Porcelain Deluxe

Technical Specification :

Deluxe quality Buchner Funnels.

Very nicely glazed and stoutly made.

Glazed inside and out, except the rim.

Inner Diameter: 110mm

Filter Paper Size: 110mm

Outer Diameter: 127mm

Cup Height: 53mm

Stem Length: 95mm

Stem Diameter: 21mm

Total Height: 200mm

Capacity: 600ml

Measurements are approximate.

Picture is a representation of buchner funnels, please refer to item description for individual product information.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003