

Product Name : Buffer Solution pH 10.0	Product Code : NLAB-CHEMICALSIND16003
Description : Buffer Solution pH 10.0	
Technical Specification : Buffer Solution pH 10.0 Appearance: Clear colourless liquid Solubility: Miscible with water pH (at 25°C): 9.96 - 10.04 Capacity:- 500 ml, 2.5 Ltr, 5 Ltr	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003