

<b>Product Name :</b> Buffer Solution pH 12.0	<b>Product Code :</b> NLAB-CHEMICALSIND16004
<b>Description :</b> Buffer Solution pH 12.0	
<b>Technical Specification :</b> Prepared with distilled water and checked for accuracy. Solutions are colorless unless stated otherwise.  Capacity: 500 ml, 2.5 Ltr, 5 Ltr  pH: 12.0  Grade: Laboratory Reagent	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003