

| | |
|--|---|
| Product Name : Buffer Solution PH 2.0 | Product Code : NLAB-CHEMICALSIND16005 |
| Description : Buffer Solution PH 2.0 | |
| Technical Specification : Buffer Solution pH 2.0 Grade: Laboratory Reagent Capacity: 500 ml, 2.5 Ltr, 5 Ltr Density: 1.003 g/mL at 20 °C pH: 2.00±0.05 | |

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003