

Product Name :
Bulk Density Kit

Product Code :
NLAB-VOCATIONALTVET17006



Description :

Bulk Density Kit

Technical Specification :

A core 57mm (2.25") in diameter is extracted and held in a brass cylinder.

The cylinder and soil sample can then be placed in a pressure plate extractor or Tempe Cell apparatus.

The water-holding characteristics of the sample can be determined.

The cylinder can be used to provide a sample of known volume, allowing the bulk density to be determined.

The sampler is supplied with two wedge coring tips, driving hammer, core extractor, spanner and strap wrench for replacing coring tips, six cylinder caps, and five brass cylinders; one 60mm long, two 30mm long, and two 10mm long

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

