

Product Name :
Button Press

Product Code :
NLAB-MACHINELABWK26011



Description :

Button Press

Technical Specification :

Button Press .

Indigenous make hand operated screw type button press desined for light pressing cutting, eyelting reveting, bending etc., having centre to back min; 1.1/2' max., 2-1/4", base to guide 2.1/4", dia of screw 3/4", hole in bed 3/8".

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003