

Product Name :
C Arm O.T. Table

Product Code :
NLABMEDICAL-M2760003

Description :

C Arm O.T. Table

Technical Specification :

The 5 section table top contains radiographic sheet of 12 mm and is equipped with assembled polyurethane 50 mm thick mattresses is easy to disassemble & cleanse.

Strong & rust free base of stainless steel is provided to make the working conditions pleasant for the doctors & nurses.

Head / Leg section are interchangeable.

Auto floor locking device.

Technical Constraints

Top Dimension: L 1890x 510 mm

Height Adjustment: 760 mm-990 mm

Trandelenburg / Reverse: 30° / 30°

Lateral Tilt: 20° / 20°

Kidney Elevator: 150 mm

Back Rest (up / Down): 80° / 25°

Leg Rest (up / Down): 15° / 90°

Head Rest (up / Down): 20° / 90°

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003