


Product Name : Cadaverous Injector	Product Code : NLABMEDICAL-M3180006
--	---



Description :

Cadaverous Injector India, Cadaverous Injector manufacturer India, suppliers, Medical Cadaverous Injector India, Hospital Cadaverous Injector India

Technical Specification :

Cadaverous Injector is used for injecting Formaldehyde solution in cadaverous at much higher speed than normal gravity process.

Unit is fully covered & mounted on a portable trolley having four castor wheels for easy movement.

Unit consists of one air compressor fitted with 1/2 HP Motor which is connected with a Stainless Steel Tank of 10liters capacity meant for storing and injecting the solution.

Tank is fitted with a safety valve, pressure gauge and rubber tubing having provision for injection.

Supplied completely with electric cord, plug and suitable to work on 220V AC 50Hz supply.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003