

Product Name :
Calibrating A Pressure Sensor

Product Code :
AR144



Description :

Calibrating A Pressure Sensor

Technical Specification :

Calibrating A Pressure Sensor
Calibration unit with dead-weight piston manometer and hand-operated spindle
Electronic pressure sensor with ceramic measuring cell, integrated amplifier and voltage output
Digital display for output signal
Additional pressure sensor as cutaway model
Set of weights
Transmission medium: hydraulic oil
Process schematic on front panel

Technical Data:

Pressure sensor
Measuring range: 0...2,5bar
Supply: 24vdc
Output signal: 0...10vdc
Piston manometer with pressure piston
Diameter: 12mm
Number of weights: 5
Pressure graduations:
0,5bar
1,0bar

1,5bar
2,0bar
2,5bar
Digital display: 4 1/2 digits
Hydraulic oil: HLP ISO 32
230V, 50Hz, 1 phase

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003