

**Product Name :**  
Calibration of Load Cells

**Product Code :**  
NLAB-ENGINEERINGLB52005



**Description :**

Calibration of Load Cells

**Technical Specification :**

**DESCRIPTION**

Here simple demonstration setup is provided to experiment with calibration procedure of load cells. Here using this setup one can calibrate different types of Load cells.

**SPECIFICATIONS**

1. Common Mechanical setup to fit various types of load cells.
  - Beam Type Load cell- 1 no.
  - 'S' Type Load cell - 1 no.
  - Compression Type Load cell - 1 no.
2. Common indicator with necessary signal conditioning circuit.

3. Common Weights: - Max 20 kg in steps of 1kg

4. Digital Indicator: -

- Least Count: - 0.1 kg

- Digit: - 3 1/2

5. Powder coated enclosure.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003