

Product Name :
Callender and Griffiths Bridge

Product Code :
NLAB-ELECTRICALAB380001



Description :

Callender and Griffiths Bridge

Technical Specification :

Three in one - Carry's Foster Bridge, Potentiometer & Callender & Griffith's Bridge.

It is an improvement over the potentiometer. The 50cm long slide wire of 1ohm resistance is calibrated into 100 equal parts. 14 coils are arranged in series in dials from each being of 1 ohm resistance terminals arranged for appropriate working of different instruments. slides on brass rod.

Accuracy : + 0.05%

Housed In Wooden Attractive Cabinet Made of Australian Pine Wood.

Patch Cords 4mm (As Required Nos.)

Optional Accessories:

DC Power Supply

Platinum Resistance Thermometer

Galvanometer (Analog/Digital)

Steam Boiler With Heating Element

Rheostat/Simple Potentiometer

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003