

Product Name :
Calorimeter for Cement and Concrete

Product Code :
NLAB-CEMETLAB130027



Description :

Calorimeter for Cement and Concrete

Technical Specification :

EN 196-8, ASTM C 186, BS 4550

Consists of a Dewar vessel, which is located in an insulation wooden box, an electrically driven stirrer with constant speed and a funnel. (Dewar vessel can be easily removed from the wooden box and replaced.)

Dimensions: 350 x 250 x 680 mm

Weight: 15 kg

Power: 230V / 50Hz

Does not include: Thermometer and stirring rod

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

