

Product Name :
Cam Mechanism

Product Code :
NLAB-MECHANICALAB240004



Description :

Cam mechanism

Technical Specification :

Cam mechanism: Applied Mechanics lab equipment: Cam mechanisms are non-uniform translatable gears that help convert uniform rotating drive motions into non-uniform drive motions with or without latching. Generally speaking, cam mechanism consists of cams, engaging members and the frame. The desired movement path is defined by the geometry of the cam and produced by pressing

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003