

Product Name :
Cantilever Beam Apparatus

Product Code :
NLAB-MECHANICALAB240005



Description :

Cantilever Beam Apparatus

Technical Specification :

Cantilever Beam Apparatus: Consists of a Stainless steel beam of L type shape with horizontal graduated horizontal arm 75cm and vertical arm of 20cm, with a spring balance of 10 kg capacity to make the horizontal leg horizontal under any load. Complete with one sliding slotted set weight set 100gm to 1 kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003