

**Product Name :**  
Cantilever Beam Apparatus

**Product Code :**  
NLAB-MECHANICALAB240005



**Description :**

Cantilever Beam Apparatus

**Technical Specification :**

Cantilever Beam Apparatus: Consists of a Stainless steel beam of L type shape with horizontal graduated horizontal arm 75cm and vertical arm of 20cm, with a spring balance of 10 kg capacity to make the horizontal leg horizontal under any load. Complete with one sliding slotted set weight set 100gm to 1 kg

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003