

Product Name :
Cavitation Demonstration

Product Code :
NLAB-MECHANICALAB230004



Description :

Cavitation Demonstration

Technical Specification :

Cavitation Demonstration: Fluid Mechanics lab equipment The advanced self-contained apparatus to demonstrate and observe the basic principles of cavitation and its implications on the performance of hydraulic machines and systems. The accessory consists of a circular Venturi-shaped test section manufactured from clear acrylic to enable visualization inside the section. As the flow of water increases the pressure at the throat falls in accordance with the Bernoulli equation

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003