

Product Name :
CBR Accessories

Product Code :
NLAB-TESTINGLABM109049



Description :

CBR Accessories

Technical Specification :

CBR Essential :

- Proving ring 10 kN capacity.
- Dial gauge 0.01 mm. x 25 mm. travel.

Spares:

- Rammer 2.6 kg & 4.89 kg.
- CBR mould complete with extension collar and perforated base plate
- Perforated swell plate with adjustable stem and lock nut
- Solid base plate

- Cutting collar
- Penetration piston assembly
- Metal tripod for dial gauge
- Circular spacer disc
- Surcharge weight (Annular) 2.5 kg
- Surcharge weight (Slotted) 2.5 kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003