

Product Name :
Ceiling Shadowless Surgical Operation Theatre Light

Product Code :
NLABMEDICAL-M2750001

Description :

Ceiling Shadowless Surgical Operation Theatre Light

Technical Specification :

Dome Head: 515 mm x 2

Halogen Bulbs: 24V x 70 W x 6 nos

Light Intensity at 1 mt: 1, 60,000 Lux

Height Adjustment: 1014 mm

Dichoric Glass Reflector: 2000 mm

Color Temperature: 4200 K \pm 300

Temp rise in field: 2°-3°C

Control Unit: SMPS

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003