

**Product Name :**  
Cement Curing Cabinet

**Product Code :**  
NLAB-CEMETLAB130035



**Description :**

Cement Curing Cabinet

**Technical Specification :**

The front doors are fitted with transparent glass. The humidity is maintained from 95% to saturation by water nebulizer. The temperature is maintained to  $20\pm 1^{\circ}\text{C}$  by a immersion heater and separated refrigerator unit . The four stainless steel racks of the internal frame can support the moulds with specimens and a large number of cement prisms. It can also be used to concrete cubes and other mortar specimens. Using a simple plastic pan is also possible to cure specimen in water.

The unit has to be connected to a suitable air compressor.

Power: 1700 W approx.

Internal dimensions (lxdxh): 1115x435x1500 mm

External dimensions: 1160x550x1900 mm

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003