

Product Name :
Cement Curing Cabinet

Product Code :
NLAB-CEMETLAB130035



Description :

Cement Curing Cabinet

Technical Specification :

The front doors are fitted with transparent glass. The humidity is maintained from 95% to saturation by water nebulizer. The temperature is maintained to $20 \pm 1^{\circ}\text{C}$ by a immersion heater and separated refrigerator unit . The four stainless steel racks of the internal frame can support the moulds with specimens and a large number of cement prisms. It can also be used to concrete cubes and other mortar specimens. Using a simple plastic pan is also possible to cure specimen in water.

The unit has to be connected to a suitable air compressor.

Power: 1700 W approx.

Internal dimensions (lxdxh): 1115x435x1500 mm

External dimensions: 1160x550x1900 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003