

Product Name :
Cement Mixer

Product Code :
NLAB-EDUCATIONTVET15049



Description :

Cement Mixer

Technical Specification :

with automatic and manual control according to standards NF P15-413, EN 196, 413, 459 and BS 3892 and ISO 679, two mixing speeds, capacity approximately 5 l

EN196 / 1 compliant steel whisk

Sand standardized EN196, 16 bags of 1350g

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003