


Product Name :
Center Lathe Machine

Product Code :
NLAB-VOCATIONALTVET20001



Description :

Center Lathe Machine

Technical Specification :

Height of Center 375 to 450 mm

Admit Between Center 850 to 10750 mm

Swing Over Bed 750 to 900 mm

Swing Over Cross Slide 450 to 600 mm

Compound Slide Travel 175 mm

Spindle Sleeve (Taper Bore) MT-5

Spindle Bore 75 to 200 mm

Taper Bore in Center Sleeve MT-5

Threads & Range (1 to 6 mm pitch) 13 Inch

Threads & Range (2 to 24 TPI) 19

Ram Travel 225 mm

Bed Length 2100 to 12000 mm

Bed Width 450 mm

No. of Speeds 8 Range (RPM) 30 to 325

Cross Slide Travel 375 to 450 mm

Length of Gap in Bed 250 mm

Swing over Gap 1000 to 1150 mm

Lead Screw Diameter & Pitch 50 mm / 2TPI Metric

Motor (1440 R.P.M. / 440V. / AC-3 Phase) 5 H.P. to 7.5 H.P.

Standard Accessories:-

Steady Rest, Follow, Chuck Plate, Spindle Sleeve, Rest, Face Plate, 2 Nos. Dead Center, Change Gear Set, Tool Post Spanner, Switch Stand.

Optional Accessories:-

Rear Tool Post, Electric Motor, True Chuck, Dog Chuck, Reversible Switch, V-Belt, Electrical Coolant Pump with Fittings, Machine Lamp,

Extra Center Height, Extra Spindle Hollow.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003