

**Product Name :**  
Centisecond Timer Digital

**Product Code :**  
NLAB-EDUSCHOOLEQ174013



**Description :**

Centisecond Timer Digital

**Technical Specification :**

Mains operated electronic timer reads 0 to 99.99 seconds. The display comprises of four 7 segment digits to read seconds and centiseconds. Two pairs of colour coded 4 mm

sockets with two push buttons are controlled with a mode switch which enables either make or break pulse to be used to operate the timer. With re-set action.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003