

Product Name :
Centrifugal Pump

Product Code :
NLAB-VOCATIONALTVET16016



Description :

Centrifugal Pump

Technical Specification :

Pressure:-

Up to 3.5 kg/cm.

Discharge:-

2000 LPH to 55000 LPH.

Pipe size:-

25 x 25 mm to 80 x 80 mm.

Power rating:-

0.50 HP to 2.00 HP.

Application:-

Industries.

Irrigation.

Water supply.

Sprinkler.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003