

**Product Name :**  
Ceramic Plate Chapati Puffer

**Product Code :**  
NLAB-VOCATIONALTVET27093

**Description :**

Ceramic Plate Chapati Puffer

**Technical Specification :**

Material : Stainless Steel

Shape : Rectangular

Automation Grade : Automatic

Frequency : 60 Hz

**Features :-**

Resistant to heat

Outstanding strength

Longer operational life

Excellent performance

Precisely designed

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003