

**Product Name :**  
Chain Sharpner

**Product Code :**  
NLAB-MACHINELABWK30035



**Description :**

Chain Sharpner

**Technical Specification :**

Chain Sharpner .

- This is an essential tool for sharpening of chains.
- It ensures uniform sharpening of all cutters.
- It is very easy to operate and time saving too.
- You will saw faster and straight with correctly sharpened chain.
- Chain sharpner is recommended for better efficiency of the chain and smooth working of the chain saw machines.
- File guide stand is also available for sharpening chain.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003