

Product Name :
Chain Sharpner

Product Code :
NLAB-MACHINELABWK30035



Description :

Chain Sharpner

Technical Specification :

Chain Sharpner .

- This is an essential tool for sharpening of chains.
- It ensures uniform sharpening of all cutters.
- It is very easy to operate and time saving too.
- You will saw faster and straight with correctly sharpened chain.
- Chain sharpner is recommended for better efficiency of the chain and smooth working of the chain saw machines.
- File guide stand is also available for sharpening chain.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003