

Product Name :
Charge Sensor

Product Code :
NLAB-VOCATIONALTVET28006



Description :

Charge Sensor

Technical Specification :

Maximum input: ± 150 V

Input capacitance: 0.01 μ F

Typical bias input current: 0.005 pA

Instrument time constant: 0.1 s

Ranges:-

± 2 V (± 20 nC)

± 0.5 V (± 5 nC)

± 10 V (± 100 nC)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003