

**Product Name :**  
Chemical Reactor Training Equipment

**Product Code :**  
NLAB-ENGINEERINGLB36008



**Description :**

Chemical Reactor Training Equipment

**Technical Specification :**

The contact time and reaction temperature can be varied.

The unit contains two supply tanks for the liquid chemicals.

The properties of different chemical reactors can be investigated by performing a specification reaction using this trainer.

Their influence on the conversion of reactants is observed in the respective reactor.

The reaction is moderated using a controlled hot water circuit.

The reactors are available as accessories.

Two variable speed diaphragm pumps pump the reactants through the reactor.

Reactant conversion and reaction temperature are measured with a combined conductivity and temperature measuring unit.

It is also possible to use PC data acquisition.

The resulting reaction product is collected in a stainless steel tray in the unit.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003