

**Product Name :**  
Chemical Reactors Bench

**Product Code :**  
NLAB-ENGINEERINGLB36029



**Description :**

Chemical Reactors Bench

**Technical Specification :**

Digital conductivity meter with the sensor.

Electrical control unit with earth leakage circuit breaker, 2 displays for r.p.m. of pumps, a display for stirring rate and an electronic thermostat.

Base Unit

2x tanks for feeding the reagents with a capacity of 5 litres each for feeding reagents

Circulation pump for hot water with electronic temperature control and safety thermostat manufactured from stainless steel.

2x peristaltic feed pumps with a flow rate of 0 to 118 millilitres per minute

We are well-known manufacturers, OEM suppliers of Chemical Reactors Bench for Chemical Engineering Lab Equipments. Contact us for high quality Chemical Reactors Bench for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003