

Product Name : Chiller Water Bath	Product Code : NLAB-EDUCATIONTVET16019
Description : Chiller Water Bath	
Technical Specification : Chiller Water Bath - Some of the other strengths of the system comprise circulation pump with 15 LPM capacity, microprocessor-based PID digital temperature Indicator-cum-controller among others. Refrigerated (Chiller) Water Bath is available with the outer body constructed using powder coated mild steel sheets. Further, the presence of inner bath chamber with lids made of stainless steel as well as PUF insulation provided between spaces of inner and outer chambers also support accurate functioning and minimum thermal loss.	
Application:-	

For analytical, diagnostic and research purposes.

For thermal control applications.

Specification:-

Cooling Rate Range	
Volume	
Below +5 °C	
At Ambient To 0°C	
At Ambient To -10°C	
At Ambient To -20°C	

Value	Value
Value	Value
Value	Value
Value	Value

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003