

Product Name : Chiller Water Bath	Product Code : NLAB-EDUCATIONTVET16019
---	--

--

Description :

Chiller Water Bath

Technical Specification :

Chiller Water Bath -

Some of the other strengths of the system comprise circulation pump with 15 LPM capacity, microprocessor-based PID digital temperature Indicator-cum-controller among others. Refrigerated (Chiller) Water Bath is available with the outer body constructed using powder coated mild steel sheets.

Further, the presence of inner bath chamber with lids made of stainless steel as well as PUF insulation provided between spaces of inner and outer chambers also support accurate functioning and minimum thermal loss.

Application:-

--

For analytical, diagnostic and research purposes.

For thermal control applications.

Specification:-

Cooling Bath Range	
Volume	
Below +5 °C	
Above 0°C	
Above -10°C	
Above -20°C	

Chamber To -30°C	
Chamber To 0°C	
Chamber To -10°C	
Chamber To -20°C	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003